sysorb®

Bioresorbierbare
Interferenzschraube

Operationstechnik

sysorb®

swiss made

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5. Publikationen

Pina S., Ferreira J.M.F. (2012)

Bioresorbable Plates and Screws for Clinical Applications: A Review. *Journal of Healthcare Engineering* 2(3): 243-260

Ettinger M., Liidakis E., Haasper C., Hurschler C., Breitmeyer D., Krettek C., Jagodzinski M. (2011)

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Luckachan G. E., Pillai C. K. S. (2011)

Biodegradable Polymers – A Review on Recent Trends and Emerging Perspectives. *Journal of Polymers and the Environment* 19(3): 637-676

Konan S., Haddad F.S. (2009)

A clinical review of bioabsorbable interference screws and their adverse effects in anterior cruciate ligament reconstruction surgery. *The Knee* 16(1): 6–13

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Hunt P., Rehm O., Weiler A. (2006)

Soft tissue graft interference fit fixation: observations on graft insertion site healing and tunnel remodeling 2 years after ACL reconstruction in sheep. *Knee Surg Sports Traumatol Arthrosc* 14(2): 1245-1251

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Thanki P. N., Dellacherie E., Six J.-L. (2005)

Prevailing mechanisms of the hydrolytic degradation of oligo(d,l-lactide)-grafted dextrans. *European Polymer Journal* 41(7): 1546–1553

Jayamoorthy T., Field J.R., Costi J.J., Martin D.K. (2004)

Biceps tenodesis: a biomechanical study of fixation methods. *Journal of Shoulder and Elbow Surgery* 13(2): 160-164

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Biomechanical Comparison of Hamstring and Patellar Tendon Graft Anterior Cruciate Ligament Reconstruction Techniques: The Impact of Fixation Level and Fixation Method Under Cyclic Loading. *Arthroscopy* 18(3): 304–315

Weiler A., Unterhauser F.N., Bail H.J., Hüning M., Haas N.P. (2002)

Smooth muscle actin is expressed by fibroblastic cells of the ovine anterior cruciate ligament and its free tendon graft during remodeling. *Journal of Orthopaedic Research* 20(2): 310–317

Weiler A., Peine R., Pashmineh-Azar A., Abel C., Südkamp N. P., Hoffmann R. F. G. (2002)

Tendon Healing in a Bone Tunnel. Part I: Biomechanical Results After Biodegradable Interference Fit Fixation in a Model of Anterior Cruciate Ligament Reconstruction in Sheep. *Arthroscopy* 18(2): 113–123

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Biomechanical Properties and Vascularity of an Anterior Cruciate Ligament Graft Can Be Predicted by Contrast-Enhanced Magnetic Resonance Imaging. A Two-Year Study in Sheep. *Am J Sports Med* 29(6): 751-761

McGuire, D. A., Barber F. A., Milchgrub, S., Wolchok J. C. (2001)

A postmortem examination of poly-L lactic acid interference screws 4 months after implantation during anterior cruciate ligament reconstruction. *Arthroscopy* 17(9): 988–992

Stähelin A. C., Südkamp N. P., Weiler A. (2001)

Anatomic Double-Bundle Posterior Cruciate Ligament Reconstruction Using Hamstring Tendons. *Arthroscopy The Journal of Arthroscopic and Related Surgery* 17(1): 88–97

John C. J., Kelly K. J., Hearn C. T., Martin D. K. (2001)

Comparison of Torsional Strengths of Bioabsorbable Screws for Anterior Cruciate Ligament Reconstruction. *Amer. J. Sports Medicine* 29: 575-580

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Brand J., Weiler A., Caborn D.N.M., Brown C.H., Johnson D. L. (2000)

Graft Fixation in Cruciate Ligament Reconstruction. *Amer. J. Sports Medicine* 28: 761-774

Weiler A., Hoffmann R.F.G., Stähelin A. C., Helling H.-J., Südkamp N. P. (2000)

Biodegradable Implants in Sports Medicine. The Biological Base. *Arthroscopy* 16(3): 305-321

Middleton J.C., Tipton A.J. (2000)

Synthetic biodegradable polymers as orthopedic devices. *Biomaterials* 21: 2335-2346

Stadelmaier D.M., Lowe, W.R., Ilahi, O.A., Noble, P.C., Kohl H.W. (1999) Cyclic pull-out strength of hamstring tendon graft fixation with soft tissue interference screws influence of screw length. *Am J Sports Med* 27(6): 778-783

Schwacha G., Vert M. (1999)

In vitro and in vivo degradation of lactic acid-based interference screws used in cruciate ligament reconstruction. *International Journal of Biological Macromolecules* 25 (1-3): 283-291

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Hoffmann R. F. G., Peine R., Bail H. J., Südkamp N. P., Weiler A. (1999)

Initial Fixation Strength of Modified Patellar Tendon Grafts for Anatomic Fixation in Anterior Cruciate Ligament Reconstruction. *Arthroscopy* 15: 392-399

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Hamstring Tendon Fixation Using Interference Screws: A Biomechanical Study in Calf Tibial Bone. *Arthroscopy* 14: 29-37

Weiler A., Windhagen H. J., Raschke M. J., Laumeyer A., Hoffmann R.F.G. (1998)

Biodegradable Interference Screw Fixation Exhibits Pull-Out Force and Stiffness Similar to Titanium Screws. *Amer. J. Sports Medicine* 26: 119-128.

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Clinical Degradation and Biocompatibility of Different Bioabsorbable Interference Screws: A Report of Six Cases. *Arthroscopy* 13: 238-244

Ignatius A.A., Claes L. E. (1996)

In vitro biocompatibility of bioresorbable polymers: poly(L, DL-lactide) and poly(L-lactide-co-glycolide). *Biomaterials* 17: 831-839

Weiler A., Hoffmann R.F.G., Stähelin A. C., Bail H.J., Siepe C. J., Südkamp N. P. (1996)

Hamstring tendon Fixation using Interference Screws: A Biomechanical Study in Calf Tibial Bone. *Biomaterials* 14: 29-37

Stähelin A.C., Liedtke H. (1995)

A Bioresorbable Screw for Arthroscopy: From Polymer to Product. *Medical Plastics* 9: 8.1-8.3

Stähelin A.C., Feinstein R., Friederich N.F. (1995)

Clinical experience using a bio-absorbable interference screw for ACL reconstruction. *Ortho Trans* 19: 287-288